

SLFP5030



This is an excellent solution for “printing underneath the silicone”. This liner is designed to back print a message in mirror image to read right through the silicone. You can then cover the mirror image with an opaque white. A varnish over the white would help printer feed and help keep the white from off setting.

SLF5030 1.2 mil Clear PET liner med. release

TC50303 Water Gloss Primer

Ship Date:	Min Order:	Master Width:	Liner Caliper:	Overall Caliper:	Shipping Weight:
2 weeks	3,000 msi	54"	1.2 mils	1.3 mils	60#/mmsi