

# SLFP5030



This is an excellent solution for “printing underneath the silicone”. This liner is designed to back print a message in mirror image to read right through the silicone. You can then cover the mirror image with an opaque white. A varnish over the white would help printer feed and help keep the white from off setting.

**SLF5030 1.2 mil Clear PET liner med. release**

**TC50303 Water Gloss Primer**

| Ship Date: | Min Order: | Master Width: | Liner Caliper: | Overall Caliper: | Shipping Weight: |
|------------|------------|---------------|----------------|------------------|------------------|
| 2 weeks    | 3,000 msi  | 54"           | 1.2 mils       | 1.3 mils         | 60#/mmsi         |