

TSG065



This a two liner transfer tape for use in making integral labels for blood bags. . Key Selling Points 1. The adhesive is FDA compliant for blood bag labels. This product should be used when our label is the primary label affixed to the blood bag. Secondary labels do not require the same premium adhesive. 2. The form stock chosen should simulate a latex face stock. 3. New orders should be subject to test. A test should include the laser printing and the application of the label to the blood bag. 4. Please note that applying the label at room temperature is best. Any condensation on the bag could affect the adhesive performance.

SL1055 42# SCK Liner - Very Easy Release

TC50619 All temp perm adhesive . Blood bag compliant. Excellent low temperature performance.

SL10040A 2.5 Mil Clay Coated Tight Release Laser Liner

Ship Date:	Min Order:	Master Width:	Liner Caliper:	Overall Caliper:	Shipping Weight:
1 Week	1,000 msi	51.5"	2.5 mils	5.5 mils	219#/mmsi